ARTech Laboratory, Inc. Checklist

Ear

Please read "Photographic Policy"
Photos taken of both sides before casting.
Please read "Ear Casting Procedure"
Lab Plaster (not plaster bandage) impression of sound side
Dental stone impression of residual affected side.
Color choices and Measurements filled out on Work Order

<u>If patient has had a prosthesis previously, it would be helpful to us to have it, or at least a photo of it.</u> They usually want the new prosthesis to look similar to the one they are using now unless there is some problem with the current one.

We do not need a lot of photos that are different colors. One or two prints of each angle, if they are color-true, will be sufficient.

It is not necessary to return the color chart with the models, unless it was printed off our web site.